

3V XTRA

Aspiration Catheter

CE₀₀₈₆

S3V
VASCULAR TECHNOLOGIES

Mfg. Lic. No.
KTK/28/413/2013

3V-XTRA Aspiration catheter provides you with the best solution for thrombus removal
3V-XTRA is a dual lumen catheter, with a distal radiopaque tip marker and proximal luer-lock port
Smooth tip designed to increase deliverability

Swift Operation

Best in class
Aspiration Speed

Optimal trackability

Unique hydrophilic Coating

Bounteous Lumen Area

Extensive Lumen area
available for extraction

Excellent Pushability

With an
innovative shaft

Kink resistance

Specifically designed for
aspiration catheter

Uncompromised Quality
Unsurpassed Value

For more information, please visit us at www.s3vvascular.com

3V XTRA

Aspiration Catheter

CE₀₀₈₆

S3V
VASCULAR TECHNOLOGIES

Mfg. Lic. No.
KTK/28/413/2013

Technical Specifications

Delivery System	Rapid Exchange
Rapid Exchange Length (cm)	25 cm
Catheter Length (cm)	145 cm
Max. Guide Wire (in.)	0.014" (OD ≤ 0.38 mm)
Guiding catheter / introducer with minimum inner lumen diameter	For 6F: 0.071" (1.80 mm)
	For 7F: 0.081" (2.06 mm)
Aspiration Time (20 mL of water)	For 6F: max 18s
	For 7F: max 10s

Guiding catheter
compatibility

3V Xtra 6F
3V Xtra 7F

Ordering Information

TYPE / MODEL	PRODUCT CODE
XTRA HYDRO 6F S Stiff	X6FS
XTRA HYDRO 7F S Stiff	X7FS

Contact in India:

S3V Vascular Technologies (P) Ltd.
#772, Balaji Chambers,
New Kalidasa Road, Vijaynagar 1st Stage,
Mysore - 570 018, Karnataka, INDIA.
Tel: +91 821 2410016
E-mail: info@s3vvascular.com

Uncompromised Quality
Unsurpassed Value

For more information, please visit us at www.s3vvascular.com